


SaeboStretch® Measurement Form

Dynamic Hand Splint



Indications

- Has mild to moderate tone or soft tissue shortening.
- Can passively achieve neutral wrist extension with composite finger extension with minimal effort.



With fingers adducted, measure straight across the proximal phalanges (2-5) proximal to the PIPs. Avoid wrapping the measuring tape.

Small: 2 1/4" - 2 3/4"
Medium: 2 3/4" - 3 1/4"
Large: 3 1/4" - 4"

Hand: Left: Right:

Size: Small: Medium: Large:

If undecided about the size,
assess finger length.
Longer fingers? Consider a larger size.

